



Selected Foods & Selected Environmental Surfaces:

The **Applied Biosystems MicroSEQ® *Listeria* spp. Detection Kit #021108** has been validated and certified as a *Performance Tested Methodsm* by the AOAC Research Institute as an effective method for the detection and identification of *Listeria* spp. in a wide variety of foods and environmental surfaces.

Dairy: Pasteurized whole cows milk, dry infant formula

Meat: Hot dogs, roast beef

Seafood: Lox (smoked salmon)

Environmental Surfaces: Stainless steel, plastic cutting board, ceramic tile, rubber sheets & concrete sealed with Seal Hard®

481 N. Frederick Ave., #500, Gaithersburg, MD 20877 USA
301-924-7077 x 129